

Circle of Willis CTA

Circle of Willis CTA – Workflow Overview


Overview:

The Circle of Willis is easily post-processed using Automated bone removal and the Vessel Grow tool. Analysis of any vessel can be accomplished using the Vessel Probe tool. This application tracks the centerline of the vessel and produces Curved Planar Reformations. The Cross Vessel view displays the lumen area and diameter of the vessel. The views can be edited, measured, captured and exported.

The Steps:

- **Load** the data and select the appropriate protocol **Vascular CoWCT**.
- **Pick** 3D Analysis.
- **Review** images using 5-on-1 screen format including the Point of Interest (POI) Cube.
- **Verify** Segmented Vessels are displayed.
- **Add** vessels using Vessel Grow.
- **Edit** the Vessel Grow Area.
- **Select** the Vessel Probe tool and click a point inside the vessel.
- **Add** additional vessels by choosing the vessel using the Vessel Probe tool.
- **Use** the Extend tool to extend the internal carotid or vertebral arteries.
- **Name** the vessels.
- **Check** for accuracy and **Edit** Centerline.
- **Select** Native thickness MPR to view a VR MIP image.
- **Create** 3D Batch Rotations of vessels.
- **Create** and **Export** snapshots and batch reformats.

Circle of Willis CTA – Select the Application




Circle of Willis CTA – Bone Removal

**Select Protocol Vascular
CoW CT.**

Select 3D Analysis.

Click on Pick




The screenshot shows the software interface for a Circle of Willis CTA. On the left, a 'Protocol' list includes various CT protocols, with 'Vascular: CoW CT' selected. Below the list is a 'Back To Viewer' button and a checked checkbox labeled 'Automatically remove bone'. On the right, a 3D viewer displays a brain scan with the Circle of Willis highlighted in orange. The viewer has a '2D' dropdown menu, a 'Pick' button in the top right, and a '3D Analysis' dropdown menu with another 'Pick' button in the bottom right. A red arrow points to the 'Pick' button in the bottom right of the 3D viewer.


Tip: Preview image from the Gallery Tab. If the contrast is not dense enough the vessels will not display well.
Clear the box next to Automatically Remove Bone to restore bone to the image.

Circle of Willis CTA – POI Cube

The POI Cube shows the volume surrounding the current crosshair.




Circle of Willis CTA – Segment Vessels




Circle of Willis CTA – Edit Vessel Grow

The steps below will help edit the blue selection area and remove unwanted anatomy.




Circle of Willis CTA – Vessel Probe Analysis

Vessel Probe is an option for CoW to isolate and analyze the vessels.




Extend will not go beyond the point you select.


Circle of Willis CTA – Rename the Vessel



Circle of Willis CTA – Centerline Edits




Circle of Willis CTA – Snapshot




Tip: It is important to remember to take a snapshot. **Only a snapshot** can be restored at a later date. Once a snapshot is restored, you will be able to continue post-processing your image.

Circle of Willis CTA – Native Thickness MPR



Tip: This tool provides the capability to view the CoW in a VR MIP projection.

Circle of Willis CTA – Native Thickness MPR





Tip: Locate the ROI. Press the wheel to display the yellow crosshairs. Center the ROI on the crosshairs. This will lock the ROI at the centerpoint and allow you to rotate the surrounding vessels.

Circle of Willis CTA – Front Cut Plane


Right mouse click in 3D image and click Front Cut Plane.

Hold down right mouse and move up/down to cut into the Volume Rendered image.


To rotate image – select the 3D tool and use left mouse to rotate.



Circle of Willis CTA – 3D Batch Rotation



Circle of Willis CTA – Export



Study List

Search: default

Filter By: Patient ID, Patient Name, Date of Birth, Accession Number, Study Description, Modality, Study Date

Status	Patient Name↑	Patient ID	Sex	Modality	# Images	Study Date/time
	CARDIAC IRO, 236, CARDIAC IRO, 236	CA1V-CUH-236	F	CT	725	17 Jul 2015 09:46
	CAROTID CTA SUBT DLP 279	AA1-NH-005	F	CT	1247	30 Jun 2015 12:47
	CCTA 0.35 mSv	CA1V-CV-106		CT	1	16 Jul 2015 05:19
	Chest - Lung Subtraction	EA1V-TUH-104	M	CT	1826	23 Mar 2015 10:25
	CHEST BREAST CA, XLFOV	ULB16-TH-026	F	CT	451	12 May 2015 09:12
	Chest c+, emphysema	EL16-TH-003	M	CT	935	24 Mar 2015 10:30
	CHEST CTA	EL16-TH-001	F	CT	384	30 Apr 2015 08:21
	CHEST LG PT, BONY LESION	UA1-TUH-096	F	CT	2733	15 Dec 2014 13:41
	COLON CA NCER, NECROTIC TUMOR	UCELB-WH-032	F	PT, CT	444	29 May 2015 07:01
	COLON CXL 102	CM64-AH-102		CT	2907	02 Jun 2015 04:13
	COLON S&P 1, DLP 293	CA1V-AH-301		CT	2903	09 Oct 2014 06:01
	COW CTA SEMAR	AA1V-HV-011	F	CT	1025	13 Jan 2015 13:06
	COW CTA SUBT 012	AA1V-HV-012	M	CT	1183	24 Mar 2015 10:39
	COW CTA SUBT RT MCA	AA1V-HV-013	F	CT	721	01 Apr 2015 07:58
	COW RCA ANEURYSM	CA1-HV-102		CT	962	27 Feb 2015 09:59
	COW sub.	EA1N-HV-010	F	CT	844	29 Nov 2014 14:06
	CSPINE #MTN BIKE	AA1-NH-004	M	CT	1420	05 May 2015 14:28
	CTA AAA	UM80-AH-042	F	CT	1492	20 Jun 2014 07:29
	CTA AAA_103	CM64-AH-103		CT	1274	10 Jul 2015 09:53
	CTA Carotids + stent	EL16-NH-001	F	CT	1319	24 Mar 2015 11:42
	FACIAL BONES_100	CM64-HH-100		CT	917	11 Jul 2015 04:31
	HEAD/NECK C&F/U, 426 LBS	UCELB-WH-044	M	CT, PT, ...	546	10 Aug 2015 09:15
	HEAD/NECK C&F/U, 426 lbs	UCELB-WH-044	M	PT, CT	492	20 May 2015 06:31
	Hematuria 403	CM80-0H-403	M	CT	1443	17 Mar 2015 07:19

Applications | **Results** | Study

COW CTA SUBT 012 | AA1V-HV-012 | CT Angiogram COW | 24 Mar 2015 10:39

CTA/SUB-BRAIN/vitreuser 24 Aug 20... 1

CTA/SUB-BRA/vitreuser 24 A

- Delete
- Publish
- Export**
- Save to Media

Circle of Willis CTA – Export

Export locations are listed in the Destination section. Choose desired destinations, click Export.

Export

Content

- 1 Studies selected
- Add results and reports

1 Studies | 0 Results | 0 Reports will be exported

Destination

- VNA
- PACS
- lungcad
- lungcad

Options

- Include annotations on images
- Convert to grayscale
- Reduce size
- Remove Vitrea results flag

Edit Patient or Anonymize

- Edit Patient
- Anonymize

Export **Cancel**

Circle of Willis CTA – Workflow Summary

Summary:

Selecting the **Vascular COW Application** Protocol you can:

- **Verify Segmented Vessels** are displayed.
- **Add** vessels using **Vessel Grow**.
- Select the **Point-of-Interest mode** to investigate surrounding vessels.
- Select the **Vessel Probe** tool to interrogate vessels.
- Select **Extend** to continue the artery down to the internal carotid arteries or vertebral arteries.
- **Name** each vessel in the Vessel Management area.
- **Edit** the Centerline of a vessel probe.
- Create **Batch Rotations** of the Circle of Willis.
- **Select** Native thickness MPR to view the VR MIP image.
- **Export** Images to PACS or other DICOM destinations from the Results Tab.

V i T A L®

Canon
CANON GROUP

Vital is a registered trademark in the U.S. and may have protection in other countries.
© 2019 Vital Images, Inc. | www.vitalimages.com